


Código	Espessura	Largura (L)	Comprimento (C)
LSM016	0.40 mm	10.00 mm	16.00 mm
LSM017	0.40 mm	8.00 mm	13.00 mm
LSM018	0.40 mm	6.00 mm	10.00 mm
LSM019	0.40 mm	6.00 mm	13.00 mm
LSM020	0.40 mm	8.00 mm	17.00 mm
LSM021	0.40 mm	10.00 mm	20.00 mm



Código	Acople	Espessura	Ø Haste	Largura (L)	Comprimento (C)	Comprimento (CC)
LRM001-A	A	0.50 mm	3.175 mm	6.00 mm	65.00 mm	15.00 mm
LRM002-A	A	0.50 mm	3.175 mm	6.00 mm	65.00 mm	34.00 mm
LRM003-A	A	0.50 mm	3.175 mm	6.00 mm	65.00 mm	22.00 mm
LRM004-A	A	0.50 mm	3.175 mm	6.00 mm	65.00 mm	27.00 mm
LRM005-A	A	0.50 mm	3.175 mm	9.00 mm	65.00 mm	16.00 mm
LRM006-A	A	0.50 mm	3.175 mm	9.00 mm	65.00 mm	25.00 mm
LRM007-A	A	0.50 mm	3.175 mm	5.00 mm	66.00 mm	17.60 mm
LRM008-A	A	0.50 mm	3.175 mm	5.00 mm	66.00 mm	22.00 mm
LRM009-A	A	0.50 mm	3.175 mm	5.00 mm	69.00 mm	24.00 mm