

BROCAS CIRÚRGICAS

Brocas Esféricas Cortantes - Média Remoção

Fresas redondas destinadas a cortes ósseos e acabamentos moderados.



Registro ANVISA 10309310020

Código	Descrição do produto	Diâmetro (ØD) mm	Comprimento total (L) mm
03092000001	BROCA ESFÉRICA MR M Ø 1,0 X 44,5 MM	1,0	44,5
03092000002	BROCA ESFÉRICA MR M Ø 1,5 X 44,5 MM	1,5	44,5
03092000004	BROCA ESFÉRICA MR M Ø 2,0 X 44,5 MM	2,0	44,5
03092000005	BROCA ESFÉRICA MR M Ø 2,5 X 44,5 MM	2,5	44,5
03092000006	BROCA ESFÉRICA MR M Ø 3,0 X 44,5 MM	3,0	44,5
03092000007	BROCA ESFÉRICA MR M Ø 3,5 X 44,5 MM	3,5	44,5
03092000008	BROCA ESFÉRICA MR M Ø 4,0 X 44,5 MM	4,0	44,5
03092000009	BROCA ESFÉRICA MR M Ø 4,5 X 44,5 MM	4,5	44,5
03092000010	BROCA ESFÉRICA MR M Ø 5,0 X 44,5 MM	5,0	44,5
03092000011	BROCA ESFÉRICA MR M Ø 5,5 X 44,5 MM	5,5	44,5
03092000012	BROCA ESFÉRICA MR M Ø 6,0 X 44,5 MM	6,0	44,5
L = 44,5 mm - Compatíveis com Ponta Drill (Cód. HM011100) e Ponta Drill Angular Multiplicadora (Cód. HM011500)			
03092000027	BROCA ESFÉRICA MR M Ø 1,0 X 70 MM	1,0	70
03092000028	BROCA ESFÉRICA MR M Ø 1,5 X 70 MM	1,5	70
03092000030	BROCA ESFÉRICA MR M Ø 2,0 X 70 MM	2,0	70
03092000031	BROCA ESFÉRICA MR M Ø 2,5 X 70 MM	2,5	70
03092000032	BROCA ESFÉRICA MR M Ø 3,0 X 70 MM	3,0	70
03092000033	BROCA ESFÉRICA MR M Ø 3,5 X 70 MM	3,5	70
03092000034	BROCA ESFÉRICA MR M Ø 4,0 X 70 MM	4,0	70
03092000035	BROCA ESFÉRICA MR M Ø 4,5 X 70 MM	4,5	70
03092000036	BROCA ESFÉRICA MR M Ø 5,0 X 70 MM	5,0	70
03092000037	BROCA ESFÉRICA MR M Ø 5,5 X 70 MM	5,5	70
03092000038	BROCA ESFÉRICA MR M Ø 6,0 X 70 MM	6,0	70
L = 70 mm - Compatíveis com Ponta Drill Média (Cód. HM011200) e Ponta Drill Multiplicadora (Cód. HM011400)			
03092000040	BROCA ESFÉRICA MR M Ø 1,0 X 115 MM	1,0	115
03092000041	BROCA ESFÉRICA MR M Ø 1,5 X 115 MM	1,5	115
03092000043	BROCA ESFÉRICA MR M Ø 2,0 X 115 MM	2,0	115
03092000044	BROCA ESFÉRICA MR M Ø 2,5 X 115 MM	2,5	115
03092000045	BROCA ESFÉRICA MR M Ø 3,0 X 115 MM	3,0	115
03092000046	BROCA ESFÉRICA MR M Ø 3,5 X 115 MM	3,5	115
03092000047	BROCA ESFÉRICA MR M Ø 4,0 X 115 MM	4,0	115
03092000048	BROCA ESFÉRICA MR M Ø 4,5 X 115 MM	4,5	115
03092000049	BROCA ESFÉRICA MR M Ø 5,0 X 115 MM	5,0	115
03092000050	BROCA ESFÉRICA MR M Ø 5,5 X 115 MM	5,5	115
03092000051	BROCA ESFÉRICA MR M Ø 6,0 X 115 MM	6,0	115
L = 115 mm - Compatíveis com Ponta Drill Longa (Cód. HM011300)			

